



PRODUCT FEATURES & DESCRIPTIONS

USER REGISTRATION

- **EHR connected invitation**
 - Pull information from the EHR
- **Non-EHR connected invitation**
 - Enter all information in the Portal
- **Self-registration**
- **Self-registration with approval**
 - Self-registration but with either email domain approval or admin approval
- **Bulk registration**
 - Admins to bulk-invite HCP users

AUTHENTICATION AND PROFILE MANAGEMENT

- Email, phone number, username as the identifiers
- Forgot password workflow
- **Manage profile**
 - User information
 - Profile photo
 - Change password
 - Set default language
- **Single sign-on**
 - From other systems, such as an EHR

MESSAGING BETWEEN CLINICIAN AND PATIENT

- **Send and receive messages to/from the patient mobile app**
 - Not for real-time chat
 - Future: Video and voice calls
- **Threaded view of messages**
 - Gmail style
- **Attachments in messages**

RECORD NOTES AND REIMBURSEMENT SESSIONS

- **Capture notes for patients that are shared with the care team**
 - Notes are automatically shared with all HCPs that have access to the patient
- **Record CPT code reimbursement times**
 - Current: HCP users will enter the time spent for reimbursement purposes manually
 - Future: Web portal will automatically track the times that are eligible for reimbursement, replacing current manual reimbursement logs.

MANAGE PATIENTS

- **At-a-glance view of patients that require attention**
 - Sorted by the urgency of patient alerts
 - Priority for high-severity alerts, followed by medium- and low-severity
- **Search, filter and sort patients**

PATIENT DETAILS

- **View and modify patient information**
 - Contact information
 - Profile information
 - Care team members
 - Alerts
 - Critical information tiles
 - Enter notes and reimbursements
 - All patient disease information data in table and chart views
 - Patient-specific rule configuration

REIMBURSEMENT MANAGEMENT

- **Track session times**
 - Support for RPM CPT codes
 - Future: CCM code support
- **Generate reimbursement reports**
 - Select a time period (default calendar month) to generate reimbursement reports
- **Determine if CPT codes are eligible**
 - Time tracked will be automatically aggregated and assigned to each CPT code



REPORT GENERATION

- **Population report**
 - For HCPs to get an overview of the patient population
- **Patient report**
 - For HCPs to generate a report with all patient data
- **Reimbursement reports**
 - For admins to generate a report to submit for reimbursement
- **View previously generated reports**
- **Print reports**
- **Custom date ranges**

RULES SETUP

- **Configure rules (to trigger flags/alerts) at an organization or patient level**
- **Maximum lookback of 30 days**
- **Severity of flags**
 - Critical
 - Major
 - Minor
- **Customize thresholds**
- **Customize rules for the number of logs or number of days**

MANAGE TEAM MEMBERS

- **Admins can manage their entire care team and other admins**
 - Invite and modify team members
 - Sort, search and filter team members
 - Deactivate team members with the ability to transfer patients to another team member
 - Change the member's organization

MANAGE ORGANIZATIONS

- **Admins can create, view and modify organizations**
 - Sort, search and filter organizations
 - View digital contracts
 - Invite organization admin users

ROLES, HIERARCHIES AND PERMISSIONS

- **Roles**
 - Out-of-the-box roles:
 - *Physician*
 - *Nurse*
 - *Educator*
 - *Administrator*
 - *Sponsor administrator*
 - *Business analyst*
 - Custom roles can be created with custom permissions and feature access
- **Organization hierarchies**
 - Up to 3 tiers of organization hierarchies
 - *Sponsor -> Health System -> Hospital/Clinic*
- **Permissions**
 - All features, access and configuration can be set either at an organization or role level

ALGORITHM HOSTING

- **Host simple algorithms (rules) or complicated ML-based algorithms for clinical decisions or for patients**

EHR INTEGRATIONS

- **Current capabilities include:**
 - Patient search
 - EHR-connected patient onboarding
 - Vitals
 - Medications
 - Encounters
 - Conditions
 - Appointments
 - US EHR systems only
- **2022 capabilities include:**
 - 2-way lab results
 - EHR-connected HCP onboarding
 - 2-way observations
 - Medical images

The Compliant Platform You Can Count On

BrightInsight provides the leading global platform for biopharma and medtech regulated digital health solutions.

brightinsight.com



The BrightInsight Disease Management Solution (DMS) is not a cleared or approved medical device. The use cases presented within this technical sheet are hypothetical use cases and do not make reference to a specific product, product claims, or product branding. The Rheumagrin interface that is illustrated is for example purposes only and does not represent a real product or product branding.

Copyright © 2022 BrightInsight, Inc. BrightInsight and the BrightInsight logo are trademarks of BrightInsight, Inc. Other trademarks are the trademarks of their owners.