



ePRO/QUALITY OF LIFE SURVEYS

- Respond to survey questions, such as:
 - Single select questions
 - Multi-select questions
 - Multi-select dropdowns
 - Numeric entries
 - URL entries
 - Feelings (Bad, Okay, Good)
 - Dates
 - Times
 - Questions and schedule of completion are driven from the Clinician Portal

EDUCATIONAL MATERIALS

- Educational materials can help patients better manage their condition or therapy via:
 - External links
 - Images
 - Documents
 - Videos

MEDICATION MANAGEMENT

- Therapies and medications may include:
 - Oral medications
 - Injectable therapies
 - Infusion therapies
- Works with care plans and set schedules
 - With customizable reminders
- Prescriptions can be pulled in from the EHR or manually entered from a database
- Guidance for usage
 - IFU, videos, tutorials, etc.
- Future: Contraindications

LOG VITAL METRICS

- For each disease condition, allow patients to manually enter vital information, such as:
 - Blood glucose
 - Blood pressure
 - Weight
 - Activity
 - Nutrition
- Apple Health and Google Fit integrations

PRODUCT FEATURES & DESCRIPTIONS

AUTHENTICATION AND PROFILE MANAGEMENT

- Registration methods
 - Accept invitations from HCP users in the web portal
 - Self-registration without HCP invitation
 - Accept invitations from HCP with EHR integration
 - Self-registration without HCP invitation with EHR integration
- Biometric authentication
 - Face ID/Fingerprint for iOS, Fingerprint/Face Unlock for Android
- Patient self-profile management
 - Profile information
 - Change password
 - Modify notifications and permissions
 - Change profile photo

INTEGRATION WITH CLINICIAN WEB PORTAL

- All data collected in the patient mobile app can be transferred to the BrightInsight cloud (and therefore the web portal)
- Send and receive messages to HCPs or the entire care team

ALGORITHM HOSTING

- Drive adherence and persistence through personalized drug dosing algorithms



CHARTS AND TRENDS

- For all observations (such as blood glucose, blood pressure, weight)
 - Customizable time period
 - More information can be displayed for each data point
 - Table view displays individual measurements
- Export or save reports of historical readings to keep track or show HCPs during an office visit
 - Can also be attached to a message to send to the clinician
- Includes summaries and aggregated information

NOTIFICATIONS AND REMINDERS

- In-app notifications and reminders can be enabled for:
 - Taking medications
 - Scheduling therapies/medications
 - New messages from HCPs
 - Upcoming appointments
 - Take ePRO/QoL surveys
 - Care plan actions
 - If no measurements have been taken in X days

SYNC CONNECTED DEVICES

- BrightInsight device SDK supports 400+ Bluetooth medical devices across 14 measurement types from 50 device vendors
- Devices include:
 - Spirometers
 - Thermometers
 - Weight scales
 - Pill reminders
 - Activity trackers
 - Blood glucose meters
 - ECGs
 - Ketone meters
 - Fetal Dopplers
 - Heart rate monitors
 - Pulse oximeters
 - Sleep trackers
 - Lactate meters

EHR INTEGRATIONS

- Current capabilities include:
 - Patient search
 - EHR-connected patient onboarding
 - Vitals
 - Medications
 - Encounters
 - Conditions
 - Appointments
 - U.S. EHR systems
- 2022 capabilities include:
 - 2-way lab results
 - EHR-connected HCP onboarding
 - 2-way observations
 - Medical images

The Compliant Platform You Can Count On

BrightInsight provides the leading global platform for biopharma and medtech regulated digital health solutions.

brightinsight.com



The BrightInsight Disease Management Solution (DMS) is not a cleared or approved medical device. The use cases presented within this technical sheet are hypothetical use cases and do not make reference to a specific product, product claims, or product branding. The Rheumagrin interface that is illustrated is for example purposes only and does not represent a real product or product branding.

Copyright © 2022 BrightInsight, Inc. BrightInsight and the BrightInsight logo are trademarks of BrightInsight, Inc. Other trademarks are the trademarks of their owners.